

McGuire Manufacturing Co., Inc.

60 Grandview Court Cheshire, Connecticut 06410 (203) 699-1801 Fax (203) 699-1813

customerservice@mcguiremfg.com mcguiremfg.com

PRODUCT SPECIFICATION

Part No.

1000WCSAN, 1151WCSAN 1151AWCSAN 1149WCSAN, 1152WCSAN

Inorganic Antimicrobial Offset Sink Strainer ADA Compliant

Job Name: Brass Slip Nut Brass Slip Nut NOTE: Fits all McGuire Sink Strainers.

MODEL NO.	DESCRIPTION	A	В
1000WCSAN	Offset Only	7"	12-3/4"
1151WCSAN	Offset with 151 Brass Strainer	7"	12-3/4"
1151AWCSAN	Offset with 151A Stainless Strainer	7"	12-3/4"
1149WCSAN	Offset with 149 Brass Flat Grid Strainer	7"	12-3/4"
1152WCSAN	Offset with 152 Brass Strainer	7"	12-3/4"



A registered trademark of Component Hardware Group, Inc.

SANIGUARD® product protection is designed to protect the product itself. It is not a disinfectant and does not protect users or others against food-borne and 1152WC illness or disease. Always use good hygiene and cleaning practices.

ASME A112.

1000WC, 1151WC, 1151AWC, 1149WC and 1152WC ASME A112.18.1/CSA B125.2



Specifications:

Adjustable offset sink strainer furnished with 17 gauge seamless brass waste arm, 17 gauge seamless brass tailpiece, cast brass slip nuts and heavy cast elbow. Sink strainer shall be in compliance with CSA or other recognized testing authority and bear both manufacturer and testing mark. All metal parts shall be treated, inside and outside, with SANIGUARD® antimicrobial compound.



This product can expose you to chemicals including lead, which is known to the State of California to cause cancer, and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.